



Revolutionizing Respiratory Care: Enhanced Health and Well-being

The 10L dual-flow oxygen concentrator can provide two patients with 93%±3% pure oxygen at the same time. It adopts oil-free compressors and PSA technology to ensure stable production of pure oxygen from the atmosphere for 365 days.

Features

- Total 10LPM output
- Dual flow outlet with two bottles
- 93%±3% oxygen purity output
- Robust design in excellent quality

Application

- Large flow rate of 10 liters oxygen production, more suitable for patients with moderate to severe hypoxia.
- Provide oxygen to patients with intermittent pneumonia and meet their high-flow oxygen requirements.
- Provide sufficient oxygen for mountain climbing enthusiasts, help prevent person from acquiring altitude sickness during hike.

Technical Specification

Oxygen Flow	0.5--10L/min in 2 outlets
Oxygen Concentration	93% ± 3%
Oxygen Output Pressure	0.4--0.6bar
Low Oxygen Alarm	Yes
Oxygen Display	Yes
Nebulizer	Yes
Noise Level	<55dB
Machine Body Size	414x337x666mm
Import Carton Size	500x415x795mm
Net Weight per Unit	26kg
Import Gross Weight per Carton	29.3kg
Watts	630 W



10LPM OXYGEN CONCENTRATOR DUAL FLOW