



Revolutionizing Respiratory Care: Enhanced Health and Well-being



5L oxygen concentrator delivers 90% to 96% high O2 purity from 0.5L to 5L. It can provide stable and high O2 concentration delivery, multifunction performance and complete alarm system. The robust compressor, ease-of-use, long lifespan make it the perfect choice for home care users.

Features

- A Touch operate multi-functional units
- Provide stable 93%±3% oxygen purity output for COPD users
- Oxygen Purity Display by Numbers (Accuracy00.0%)
- With timer and child lock
- Large 6' LED Light Display all functional data
- Six security warning systems
- Accumulated total running hours

Application

The applicable people are chronic obstructive pulmonary disease (COPD) patients, hypoxic users and people who need daily oxygen therapy at home. It can also be widely used in hospitals, clinics, health centers, pet hospital, beauty salon and fitness center etc.

Technical Specification

Oxygen Flow	0.5--5L/min
Oxygen Concentration	93%±3%
Oxygen Output Pressure	0.4--0.7bar
Low Oxygen Alarm	Yes
Oxygen Display	Yes
Nebulizer	Yes
Noise Level	<43dB
Watts	320 W
Machine Body Size	380x270x580mm
Import Carton Size	460x350x680mm
Net Weight per Unit	16.5kg
Import Gross Weight per Carton	18.7kg



5L OXYGEN CONCENTRATOR

AMU UK LIMITED
M17 1SR Manchester, UK
E: info@amu-equipment.co.uk