



# First Aid at Your Side!

## DEFIBRILLATOR SPECIFICATION

### Category

- Mode
- Waveform
- Energy
- Voice Prompt
- Visual Indicator
- Charge time (new, at 25 C)
  
- Analysis Time
- Controls
- Output disabled when the patient Impedance is outside limits

### Specification

- Adult
- 3.5 inch
- 320 (RGB) \*240
- Biphasic truncated exponential
- 200-joulex max
- Multiple Languages to select
- Color LCD and LED indicators
- Less than 7 sec, to 150j Less than 3 sec to 75j
- Less than 10 sec to 200j Less than 3 sec to 50j
- 9 sec
- Three buttons - ON/OFF, Shock, Set
  
- 20Ω to 200Ω



# Automated External Defibrillator Trainer HBO-09T

**AMU UK LIMITED**  
M17 1SR Manchester, UK

E: [info@amu-equipment.co.uk](mailto:info@amu-equipment.co.uk)