



# First Aid at Your Side!


## DEFIBRILLATOR SPECIFICATION

### Category

Mode  
Waveform  
Energy  
Voice Prompt  
Visual Indicator  
Charge time (new, at 25 C)  
  
Analysis Time  
Controls  
Output disabled when the patient  
Impedance is outside limits

### Specification

Adult  
3.5 inch  
320 (RGB) \*240  
Biphasic truncated exponential  
200-joules max  
Multiple Languages to select  
Color LCD and LED indicators  
Less than 7 sec, to 150j Less than 3 sec to 75j  
Less than 10 sec to 200j Less than 3 sec to 50j  
9 sec  
Three buttons - ON/OFF, Shock, Set  
  
20Ω to 200Ω



PC Integration and Viewer Software

# Automated External Defibrillator HBO-09

**AMU UK LIMITED**  
M17 1SR Manchester, UK

E: [info@amu-equipment.co.uk](mailto:info@amu-equipment.co.uk)



# AUTOMATED EXTERNAL DEFIBRILLATOR

## HBO-09

### ENVIRONMENT SPECIFICATION

#### Category

Operating temperature  
Operating humidity  
Atmospheric Pressure  
Storage and transport temperature  
Particle and water ingress  
Drop and Shock Abuse Tolerance

#### Specification

0 °C to 50 °C  
Relative humidity between 15% and 95 %  
500hpa ~1060hPa  
-20 °C to 55 °C  
IP-55 (Battery installed)  
IEC60601 -1 -12  
(1 meter, Onto any edge, corner, Or surface while in standby mode)

### ACCESSORIES

#### Category

Battery type  
AED pads

#### Specification

Non-rechargeable, LimnO2 12V, 2.8Ah  
Disposable Adult pads



**AMU UK LIMITED**  
M17 1SR Manchester, UK

E: [info@amu-equipment.co.uk](mailto:info@amu-equipment.co.uk)